# KORONA<sup>®</sup>

## Daniela

- (D)** **DIAGNOSE-WAAGE** Gebrauchsanweisung
- (GB)** **DIAGNOSTIC SCALE** Instructions for use
- (F)** **PÈSE-PERSONNE IMPÉDANCEMÈTRE** Mode d'emploi
- (E)** **BÁSCULA PARA DIAGNÓSTICO** Instrucciones para el uso
- (I)** **BILANCIA DIAGNOSTICA** Istruzioni per l'uso
- (GR)** **ΔΙΑΓΝΩΣΤΙΚΗ ΖΥΓΑΡΙΑ** Οδηγίες χρήσης
- (RUS)** **ДИАГНОСТИЧЕСКИЕ ВЕСЫ** Инструкция по применению
- (PL)** **WAGA DIAGNOSTYCZNA** Instrukcja obsługi
- (NL)** **DIAGNOSEWEEGSCHAAL** Gebruiksaanwijzing

 Hans Dinslage GmbH  
88524 Uttenweiler, GERMANY  
Tel.: 07374/915766  
[www.korona-gmbh.de](http://www.korona-gmbh.de)

## **Inhalt**

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## Contents

1	Getting to know your instrument.....	16	8	Operation.....	22
2	Signs and symbols .....	17	9	Evaluation of results .....	24
3	Safety notes .....	17	10	Replacing batteries.....	25
4	Unit description .....	18	11	Storage and maintenance .....	26
5	Initial use .....	18	12	Disposal.....	26
6	Information .....	19	13	What if there are problems? .....	27
7	Settings .....	20			

## Included in delivery

- Diagnostic scale
- 2 x 3 V CR2032
- These instructions for use

## Dear Customer

You have chosen a quality product from Korona. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

This diagnostic scale is designed to make a significant contribution to your health.

With kind regards  
Your Korona team

## 1 Getting to know your instrument

### Function of the unit

This digital diagnostic scale is used to weigh you and to provide diagnostic data about your personal fitness.

The scale is intended for individual, household use.

The scale has the following diagnostic functions that can be used by up to 10 persons:

- Body weight measurement
- Measurement of body fat percentage
- Body water percentage
- Body mass index and
- Daily calorie requirement

The memory function lets up to 10 people access their last measured weight and the body values that are simultaneously determined.

This scale also has the following functions:

- Switch between kilograms [kg], pounds [lb] and stones [st:lb]
- Athlete mode
- Automatic shutoff function
- Battery change indicator for weak batteries

## 2 Signs and symbols

The following symbols appear in these instructions.

**Warning** Warning instruction indicating a risk of injury or damage to health.

**Important** Safety note indicating possible damage to the unit/accessory.

**Note** Note on important information.

## 3 Safety notes

### Safety notes

**Warning:**

- The scales must not be used by persons with medical implants (e.g. heart pace-makers) Otherwise their function could be impaired.
- Do not use during pregnancy.
- Do not stand on the outer edge to one side of the scale otherwise it may tip!
- Keep children away from packaging materials (risk of suffocation).

### Notes on handling batteries

**Warning:**

- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a child swallow a battery, seek medical assistance immediately.
- Batteries must not be recharged, taken apart, thrown into an open fire or short circuited.
- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable legal regulations. Do not dispose of batteries with the normal household waste.

**Warning:**

- Replace weak batteries before they discharge completely.
- Always replace all batteries at the same time and use batteries of the same type.
- Leaking batteries may damage the device. If you do not intend to use the device for longer periods, remove the batteries from the battery compartment.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.

## General notes


- The unit is for personal use only and is not intended for medical or commercial applications.
- Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
- The scales' maximum load is 150 kg (330 lb/24 st)
- Before submitting any complaint, first check the batteries and replace them if necessary.
- Repairs may only be performed by the authorised customer services team.
- All scales comply with EC Directive 2004/108/EC. Should you have any questions concerning the use of our devices, please contact your dealer or customer service.

## 4 Unit description

### Overview

- 1 Electrodes
- 2  Up button
- 3  On/Set button
- 4  Down button
- 5 Display

## 5 Initial use

### Inserting the batteries

Please remove the battery insulating strips and/or position the batteries correctly with regards to the polarity displayed in the battery compartment (as symbol  $\oplus$ ). If the scale fails to operate, remove the batteries completely and reinsert them.

### Changing the weight unit

The scale shows the values in kg. You can set the device to measure in kilograms [kg], or in pounds [lb] and stones [st:lb], using the button at the back of the device.

### Setting up the scale

Set up the scale on a level and solid surface. A solid surface is a prerequisite for a correct measurement.

## 6 Information

### The measuring principle

This scale operates according to the B.I.A. principle (bioelectric impedance analysis). This enables the measurement of physical relationships within seconds by means of an undetectable, completely harmless electric current. The body fat percentage and other physical relationships in the body can be determined by measuring the electrical resistance (impedance) and calculating constants and individual parameters such as age, height, gender and degree of activity.

Muscle tissue and water have good electrical conductivity, and therefore lower resistance. Bones and fatty tissue, on the other hand, have low conductivity as fat cells and bones conduct hardly any current as a result of their very high resistance.

Please be aware that values obtained from the diagnostic scale represent only an approximation of actual analytical medical data. Only a specialist physician can accurately determine body fat, body water, muscle percentage and bone structure using medical procedures (e.g. computed tomography).

### General tips

- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Important for the measurement: only measure body fat when barefoot and with the soles of the feet slightly moist. Completely dry soles can result in unsatisfactory measurements due to inadequate conductivity.
- Stand still and upright during the measurement.
- Wait a few hours after unaccustomed physical activity.
- Wait approx. 15 minutes after getting up in the morning to allow the water in the body to distribute.
- Remember that only long-term trends are important. Short-term changes in weight within a few days are usually caused by a loss of fluid. Body water plays an important role towards our general wellbeing.

### Limitations

When measuring body fat and other values, deviating and implausible results may occur in:

- children under approx. 10 years,
- competitive athletes and body builders,
- pregnant women (due to the amniotic fluid),
- persons with fever, undergoing dialysis treatment or with symptoms of edema or osteoporosis,
- persons taking cardiovascular medication (affecting the heart and vascular system),
- persons taking vascular dilating or vascular constricting medication,
- persons with considerable anatomic deviations of the legs in relation to total body size (length of the legs considerably shortened or lengthened).

## 7 Settings

### Setting user data

You must enter personal user parameters in order to measure your body fat percentage and other physical data.

The scale has 10 user memory locations which allows, for example, you and your family members to save personal settings and retrieve them later:

Sequence	User data	Parameters
❶	Memory position	1 to 10
❷	Sex	male  female 
❸	Age	10 to 100 years
❹	Body height	100 to 240 cm (3'-03" bis 7'-10")
❺	Athlete mode	Light/Irregular sport  Regular recreational sport  Competitive sport 


### Athlete mode

Physical exercise leads to a specific metabolism and high muscle density. This leads to an increase in the body's resistance. So that this increased resistance is not wrongly interpreted as fat, you can adapt the measurement parameters to suit different (sport-related) physical conditions. Select the level of activity that best fits your normal daily routine:

- Light/irregular sport : walking, light sport without regular training.
- Regular recreational sport : Average physical activity, approx. 30 mins 2-5 times per week (e.g. long walks, regular sport, etc.)
- Competitive sport : Intensive physical exertion, intensive training or hard physical work, at least 1 hour each day.

Always adjust the settings in all areas as follows:

- Press the buttons  and , to change a display or a value.
- Press the button , to save the setting and go to the next range of settings, or to close the settings menu.

## 8 Operation


### Weighing and measuring other body values

To simultaneously weigh yourself and determine other body values:

- Press the button  and wait until the display shows “0.0 kg”.
- Press the buttons  or , to select the desired storage space with your saved settings.

Confirm your selection with the button .

Your settings and saved values are displayed in an automatic display sequence. As soon as the display sequence has finished and “0.0 kg” is displayed again, you can step onto the scales to be weighed.

## Display of results


The results are automatically displayed in the following sequence:

		
1. Current weight	2. Calculation procedure	3. Body fat percentage (FAT %)
		
4. Body water percentage (TBW %)	5. Body mass index (BMI)	6. Calorie requirement (KCAL) in kilo calories

After displaying all values, the scales automatically switch off. You will find information on analysis on page 24.

### Weighing only

If you only want to calculate your body weight, just switch on the scales with button . Wait until "0.0 kg" is displayed.

After displaying the weight, the scales automatically switch off.

## 9 Evaluation of results

### Body fat percentage

The body fat values are shown in the display. The following body fat percentages are guide values (contact your physician for further information). The body fat chart that is also shown on the display (at the bottom) is only to provide more simple orientation.

Man						Woman					
Age	Severely under-weight	Under-weight	Normal	Over-weight	Obese	Age	Severely under-weight	Under-weight	Normal	Over-weight	Obese
≤ 30	4-8	8.1-14	14.1-20	20.1-27	27.1-70	≤ 30	4-12	12.1-17	17.1-24	24.1-31	31.1-70
> 30	4-11	11.1-17	17.1-23	23.1-30	30.1-70	> 30	4-15	15.1-20	20.1-27	27.1-34	34.1-70

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

### Body water percentage

The body water percentage is normally within the following ranges:

Man				Woman			
Age	poor	good	very good	Age	poor	good	very good
10-100	<50	50-65	>65	10-100	<45	45-60	>60

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.

### Body mass index (BMI)

The body mass index (BMI) is a number that is often called upon to evaluate body weight. The number is calculated from body weight and height, the formula is: Body mass index = body weight : height<sup>2</sup>. The unit of measure for BMI is [kg/m<sup>2</sup>]. Using the BMI, weight is classified for adults (20 years and over) with the following values (source: WHO):

Category		BMI
Underweight	Severely underweight	< 16
	Underweight	16-16.9
	Slightly underweight	17-18.4
Normal weight		18.5-25
Overweight	Overweight	25.1-29.9
Obese (overweight)	Class I obesity	30-34.9
	Class II obesity	35-39.9
	Class III obesity	≥ 40

### Daily calorie requirement

The displayed calorie requirement is the amount of energy that the active body needs each day. A person's calorie requirement increases with more physical activity and is calculated by the diagnostic scales (in athlete mode) using the selected level of activity (1-3).

To maintain the current weight, the energy that the body uses must be replaced in the form of food and drink. If less energy is taken in that is used over a long period of time, the body takes the difference from the fat stores and loses weight. On the other hand, if more energy is taken in than the calculated calorie requirement, the body cannot burn off the excess energy. The excess is stored in the body as fat and the body gains weight.

### Results in relation to time

---

Remember that only long-term trends are important. Short-term fluctuations in weight over a few days are usually the result of a loss of fluid.

---

The interpretation of the results reflects changes in total body weight, percentage of body fat and body water as well as muscle density, and takes into account the length of time over which these changes take place. Sudden changes within days must be distinguished from medium term changes (over weeks) and long term changes (months).

A basic rule is that short term changes in weight almost exclusively represent changes in water content, whereas medium and long term changes may also involve the fat and muscle percentages.

- If your weight reduces over the short term, but your body fat percentage increases or remains the same, you have merely lost water – e.g. after a training session, sauna session or a diet restricted only to rapid weight loss.
- If your weight increases over the medium term and the body fat percentage falls or stays the same, then you could have built up valuable muscle mass.
- If your weight and body fat percentage fall simultaneously then your diet is working – you are losing fat mass.
- Ideally you should support your diet with physical activity, fitness or power training. This way you can increase your muscle percentage over the medium term.
- Body fat, body water or muscle percentages should not be added (certain elements of muscle tissue also contain body water).

## 10 Replacing batteries

The batteries in the scales must be replaced if the message “LO” appears in the display. If the batteries are too flat, then the scales will switch off automatically.

---

**Note:**

- When changing batteries, always use batteries of the same type, brand and capacity.
  - Always replace all batteries simultaneously.
  - Do not use rechargeable batteries.
  - Use batteries that are free of heavy metals.
-

## 11 Storage and maintenance

The precision of the measurements and service life of the device depends on its careful handling:

### Important:


- The unit should be cleaned occasionally. Do not use abrasive detergents and never immerse the unit in water.
- Make sure that no liquids come into contact with the scales. Never submerge the scales in water. Never rinse them in running water.
- Do not stand any objects on the scales when not in use.
- Protect the scales from knocks, damp, dust, chemicals, drastic changes in temperature, electromagnetic fields and nearby sources of heat (ovens, heaters). Do not press the buttons violently or with pointed objects.
- Do not wash the unit in a dishwasher!

## 12 Disposal

Standard and rechargeable batteries should not be disposed of separately from the household waste. As a consumer, you are legally obliged to return used batteries for proper disposal. You can hand in your used batteries at public collection points in your district or sales outlets where batteries of this type are sold.

### Note:

The codes below are printed on batteries containing harmful substances:

- Pb = Battery contains lead,
- Cd = Battery contains cadmium,
- Hg = Battery contains mercury.

For environmental reasons, do not dispose of the scale with the batteries in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point.

Observe the local regulations for material disposal.

Dispose of the unit in accordance with EC Directive **2002/96/EC – WEEE** (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.

## 13 What if there are problems?

If the scale detects an error when measuring, the following is displayed.

Display	Cause	Remedy
Err	The maximum load capacity of 150 kg was exceeded.	Do not load the unit with over 150 kg.
LO	The batteries are almost empty.	Replace the batteries
- - - -	The electrical resistance between the electrodes and the sole of your foot is too high (e.g. if you have tough calluses).	Repeat the measurement barefoot. If necessary, moisten the soles of your feet slightly. If necessary, remove the calluses from the soles of your feet.
	The fat percentage is outside the measurable range (less than 4% or greater than 55%).	Repeat the measurement barefoot or, if necessary, moisten the soles of your feet slightly.

### No measurement possible

Possible errors	Remedy
The scale was not activated before you stood on it. If you step on the scale before "0.0 KG" is displayed, the scale will not work correctly.	Activate the scale correctly (wait for "0.0 KG") and repeat the measurement.
The batteries in the scale are empty.	Replace the batteries.

### Incorrect weight measurement

Possible errors	Remedy
Scale is on carpeting.	Place the scale on a level and solid surface.
Incorrect zero point of scale.	Activate the scale, wait approximately 4 seconds, repeat measurement.

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